

Made by:  
Jsm Micro inc. Montreal Canada

N°	Event	System_ID	Place	Call_Type	Prio	Zone_ID	Duration	Escalats	Escal_A
2381	0006	Cardiology	109-1	BED CALL	001	WARD-1	0000	01	0
2380	0005	Cardiology	133	BED CALL	001	WARD-1	0000	01	0
2379	0004	Cardiology	103	BED CALL	001	WARD-1	0000	01	0
2378	0003	Cardiology	102	BED CALL	001	WARD-1	0000	01	0
2377	0002	Cardiology	101	BED CALL	001	WARD-1	0000	01	0
2376	0001	CPU	---	WSX1 RESET (r...	255	---	0000	00	00
2048	1D01	3 W. Surgical	3559	TOILETTE	03	00	0000	00	00
2047	1D25	4 W. Medical	4512-1	LIT	05	00	0002	00	00
2046	1D24	4 W. Medical	4512-1	LIT	05	00	0001	00	00
2045	1D23	4 W. Medical	4524	LIT	05	00	0005	00	00
2044	1D22	4 W. Medical	4514	CORDON	02	00	0000	00	00
2043	1D21	4 W. Medical	4512-1	LIT	05	00	0001	00	00
2042	1D20	4 W. Medical	4516-2	TOILETTE	03	00	0000	00	00
2041	1D1F	4 W. Medical	4516-2	TOILETTE	03	00	0000	00	00
2040	1D1E	4 W. Medical	4516-2	TOILETTE	03	00	0000	00	00
2039	1D1D	4 W. Medical	4516-2	TOILETTE	03	00	0000	00	00
2038	1D1C	4 W. Medical	4516-2	TOILETTE	03	00	0002	00	00
2037	1D1B	3 W. Surgical	3524	TOILETTE	03	00	0000	00	00
2036	1D1A	4 W. Medical	4546	LIT	05	00	0000	00	00
2035	1D19	4 W. Medical	4546	LIT	05	00	0001	00	00
2034	1D18	3 W. Surgical	3516-1	LIT	05	00	0001	00	00
2033	1D17	3 W. Surgical	3508-2	LIT	05	00	0004	00	00
2032	1D15	4 W. Medical	4526-2	LIT	05	00	0001	00	00
2031	1D16	3 W. Surgical	3574-2	LIT	05	00	0000	00	00
2030	1D14	3 W. Surgical	3558-1	LIT	05	00	0000	00	00
2029	1D13	3 W. Surgical	3512-2	LIT	05	00	0000	00	00
2028	1D12	4 W. Medical	4544	LIT	05	00	0002	00	00
2027	1D11	4 W. Medical	4558-2	LIT	05	00	0002	00	00
2026	1D10	4 W. Medical	4558-2	LIT	05	00	0002	00	00

**Overview**

Jsm Report is a software for downloading, storing, managing and querying report information from WSX1 units. Reports are generated by the software included in the WSX1 unit .Jsm Report is customizable and can manage different types of reports depending on the installed system such as the WSX1 CTD3.0 software. Jsm Report can operate in single user mode or multi-user mode (TCP/IP network). It supports SQLite databases as well as Access databases. It can also generate automated reports.

**Features**

- manage database files containing downloaded report information from WSX1
- view data in a clear and flexible way
- quickly sort data by simply clicking on the appropriate column header
- make quick searches to quickly find something
- launch advanced queries to look for certain patterns for analysis and reporting purposes
- quickly generate print-ready reports for anormal situations, statistics, room call listing, etc.
- download report information from Jsm Micro inc. WSX1 units
- copy and export data to external programs such as Microsoft Excel
- act as a server in a network; multiple users can connect to the same database and query it
- auto-launch when Windows starts and periodically check for updated report information
- use SQLite and Access databases
- special databases and reports for logging pager calls
- generate automated reports

**Detailed call listing between 01/01/2009 and 31/12/2009**

Report generated by Julien Gascon-Samson (Jsm Micro inc.) on August 26, 2009 11:17:28 AM

**Report settings :**  
 Place : 3559  
 System ID : %  
 Zone ID : %  
 Start date : 01/01/2009  
 End date : 31/12/2009  
 Maximum displayed : 200

**Summary**  
 Total number of calls made between 01/01/2009 and 31/12/2009 for selected locations : 1 calls

**Detailed view**

N°	Place	Zone_ID	System_ID	Call_Type	Date_Time	Duration	Escalats
2048	3559	00	3 W. Surgical	TOILETTE	27-08-2009 11:46:48 AM	0000	00

(Warning : only the first 200 calls have been shown.)

**Use with the following products**

- WSX1 : CPU
- CTD/other : call management software

**Technical specifications**

- Systems : Windows XP, Vista, 7, 8, 8.1
- Platform : Microsoft .NET Framework 3.5
- Network : TCP/IP client/server, port 1003
- Security : USB Dongle or serial number
- Database : SQLite and Microsoft Access
- Connection : serial (RS232) connexion to WSX1